

## Related topics

Law of induction, magnetic flux, Maxwell's equations.

## Principle

A permanent magnet falls with different velocities through a coil. The change in the magnetic flux  $\Phi$  generates an induced voltage impulse  $U_{SS}$  is recorded with a computer interface system. Depending on the polarity of the permanent magnet the induced voltage impulse is negative or positive.

## Equipment

Cobra3 Basic Unit	12150.00	1
Power supply, 12 V <sub>-</sub>	12151.99	1
RS232 data cable, 9 poles	14602.00	1
Cobra3 Universal Plotter Software	14504.61	1
Forked light barrier, compact	11207.20	1
Support rod, round, $l = 600$ mm	02037.00	1
Boss head	02043.00	3
Tripod, -PASS-	02002.55	1
Universal clamp	37715.00	1
Glass tube, $l = 300$ mm	45126.01	1
Coil holder	06528.00	1
Coil, 600 turns, short	06522.01	1
Magnet, $d = 8$ mm, $l = 60$ mm	06317.00	1
Connecting cord, 32 A, $l = 50$ cm, red	07361.01	2
Connecting cord, 32 A, $l = 50$ cm, blue	07361.04	2
Connecting cord, 32 A, $l = 50$ cm, yellow	07361.02	1
PC, Windows® 95 or higher		





Fig. 2. Circuit diagram

Fig. 1. Experimental set-up



### Tasks

1. Measurement of the induced voltage impulse  $U_{SS}$  and the velocity of the falling magnet.
2. Evaluation of the induced voltage impulse  $U_{SS}$  as a function of the velocity of the magnet.
3. Calculation of the magnetic flux induced by the falling magnet as a function of the velocity of the magnet.

### Set-up

- According to Figs. 1 and 2.
- Adjust the light barrier in a such a manner that it is located directly above the coil. The magnet must interrupt the light beam during its fall.
- The glass tube provides the magnet with reliable guidance and also insures a stable spinless fall even at larger fall heights.

### Procedure

- Initiate recording of the measured values using parameters in accordance with Fig. 3.
- Depress the "Start measurement" button.
- Introduce the tip of the magnet in the glass tube, let it fall, catch it with one hand under the coil and manually terminate the measurement recording procedure.
- Measure the interruption period of the signal from channel IN 1 (light barrier) using the "Survey" icon (Fig. 4). In this case:  $\Delta x = \Delta t = 26.2$  ms. With a magnet length of  $l = 6$  cm, it follows that the velocity immediately before the magnet enters the coil is  $v = l / \Delta t = 2.29$  m/s.

Fig. 3. Measurement parameters





Fig. 4. Determination of the magnet's falling velocity



- Also using the "Survey" icon, measure the total amplitude  $U_{SS}$  (peak to peak) of the induction voltage (Fig. 5); in this case:  $U_{SS} = \Delta y = 2.764$  V.
- Mark unnecessary parts of the curve before and after the induction pulse and cut them out with the "Cut" icon.
- Mark the positive ( $F_1$ ) and the negative ( $F_2$ ) parts of the curve separately (cf. Fig. 7). Calculate the areas with the "Show integral" icon. Finally, add the values of the two partial areas (Fig. 6).

### Results

- For induced voltage  $U$  the following is true:

$$U = -n \frac{d\Phi}{dt}$$

where

$n$  = Number of windings in the coil

$\Phi$  = Magnetic flux

$t$  = Time

Integration provides the following:

$$\Phi = -\frac{1}{n} \int U dt = B \cdot A$$

where

$B$  = Value of the magnetic flux density

$A$  = Cross sectional area of the coil

Fig. 5. Determination of the induction voltage  $U_{SS}$





Fig. 6. Calculation of the magnetic induction



- Plot  $U_{ss}$  versus the falling velocity  $v$  of the magnet (Fig. 8). It is apparent that the induced voltage is proportional to the velocity:  $U_{ind} \sim d\Phi/dt$ .
- It follows from Fig. 9 that the integral of the induced voltage is a constant which is not a function of the velocity of field change:

$$\Phi = B A$$

## Remarks

- The curves in Figs. 4 to 6 are clearly asymmetrical, in particular for small falling heights. This is due to the non-constant velocity of the magnet on falling through the coil.
- When the polarity of the magnet reverses, the sign of the induction voltage changes and thus the course of the curve.

Fig. 8. Amplitude of the induced voltage  $U_{ss}$  as a function of the magnet's falling velocity  $v$




Fig. 7. Principle of the calculation of the magnetic induction




Fig. 9. Magnetic flux  $\Phi$ , determined from measurement curves according to Fig. 6, for differing magnet velocities

